

OXYGEN SOLUTIONS



Oxygen Concentrators,
Cylinders and more...

LF-H-10A

10 LPM High-pressure Oxygen Concentrator

- Triple filtration system ensures purity of the oxygen supply
- 24 hours non-stop working capacity
- Works best with ventilators and anesthesia machines



Suitable to use with



Ventilator



Anaesthesia machine



Hyperbaric oxygen chamber




Oxygen Inhaled Pressure Regulator

10 LPM High-pressure Oxygen Concentrator

10 LPM high-pressure oxygen concentrator is capable of delivering oxygen with a full boost in hospital emergencies and other departments. This machine supports oxygen flow for 5 patients at a time without a break. Therefore, there is no longer a need to buy individual oxygen concentrators for the patients admitted to the medical wards.



Product Features



Specifications

Name	Medical Hyperbaric Oxygen Concentrator
Model	LF-H-10A
Machine Type	Class II, Type B
Flow Rate	10L/min
Oxygen concentration	≥93% ±3%
Noise	≤55dB (A)
Voltage / Frequency	220V/50Hz, 110V/50Hz, 230V/50Hz
Rated power	1200W
Output pressure	0.32 - 0.42 Mpa
Operating Temperature	5°C - 40°C
Machine Size	640x340x656 mm
Package Size	820x440x870 mm
Net Weight	75 Kg
Gross weight	90 Kg
Working time	24 hour uninterrupted operation
Optional configuration	Diversion panels
Use with	Ventilator, anaesthesia machine, hyperbaric oxygen chambers, small oxygen filling stations



CE 0197






SANAI HEALTH


3A Sanai Health Technology Co., Ltd.

 400 187 6068

 +86 755 82727906

 www.sanaihealth.com

 3502H, Block A, Jingji Binhe Times Square, 9289
Binhe Road, Futian District, Shenzhen, China PRC.



LOVE NATURE
LOVE SOCIETY
LOVE SOCIETY