

OXYGEN SOLUTIONS



Oxygen Concentrators,
Cylinders and more...

LFY-I-5A-01 (10L Edition)

10 LPM Oxygen Concentrator

- 10 LPM with Dual Flow
- 93% medical oxygen concentration dual flow




10 LPM Oxygen Concentrator

Our home oxygen concentrator delivers all the qualities a home-care provider or a patient must be looking for. With the integrated PSA technology, this 10 LPM oxygen concentrator by 3A Sanai Health provides stable and high oxygen concentration 24 hours non-stop. Breathing purified and filtered oxygen at home is no more a problem for patients facing respiratory deficits.



Product Features



Specifications

Name	10L Medical Oxygen Concentrator
Model	LFY-I-5A-01 10L Edition
Machine type	Class II, Type B
Flow rate	10L/min
Oxygen concentration	≥93% ±3%
Noise	≤55dB (A)
Voltage / Frequency	220V/50Hz, 110V/50Hz, 230V/50Hz
Power rating	580W
Output Pressure	0.03 - 0.07 Mpa
Operating Temperature	5°C - 40°C
Machine Size	410x410x620 mm
Package Size	500x480x745 mm
Net weight	28.5 Kg
Gross weight	31.5 Kg
Working time	24 hour uninterrupted operation
Optional configuration	Nebulizer function



CE 0197





SANAI HEALTH


3A Sanai Health Technology Co., Ltd.

 400 187 6068

 +86 755 82727906

 www.sanaihealth.com

 3502H, Block A, Jingji Binhe Times Square, 9289
Binhe Road, Futian District, Shenzhen, China PRC.



LOVE NATURE
LOVE SOCIETY
LOVE SOCIETY