

OXYGEN SOLUTIONS



Oxygen Concentrators,
Cylinders and more...

LFY-I-5A-01

5 LPM Medical Oxygen Concentrators

- Equipped with 2 flowmeters and dual outlets
- Up to 93% oxygen concentration




5 LPM Oxygen Concentrator

Our oxygen concentrators are designed to filter the air with a 3-grade filtration system for 2 patients at a time. It makes sure that each one of the patients receives purified and healthy oxygen without microbial contamination. Additionally, this 5 LPM model can deliver up to 93% of purified oxygen 24 hours non-stop.



Product Features



Specifications

| | |
|-------------------------------|---------------------------------|
| Name | 5L Medical Oxygen Concentrator |
| Model | LFY-I-5A-01 |
| Machine type | Class II, Type B |
| Flow rate | 5L/min |
| Oxygen concentration | ≥93% ±3% |
| Noise | ≤50dB (A) |
| Voltage / Frequency | 220V/50Hz, 110V/50Hz, 230V/50Hz |
| Power rating | 450W |
| Output Pressure | 0.03 - 0.07 Mpa |
| Operating Temperature | 5°C - 40°C |
| Machine Size | 410x410x620 mm |
| Package Size | 500x480x745 mm |
| Net weight | 28.5 Kg |
| Gross weight | 31.5 Kg |
| Working time | 24 hour uninterrupted operation |
| Optional configuration | Nebulizer function |



CE 0197






SANAI HEALTH


3A Sanai Health Technology Co., Ltd.

 400 187 6068

 +86 755 82727906

 www.sanaihealth.com

 3502H, Block A, Jingji Binhe Times Square, 9289
Binhe Road, Futian District, Shenzhen, China PRC.



LOVE NATURE
LOVE SOCIETY
LOVE SOCIETY