

OXYGEN SOLUTIONS



Oxygen Concentrators,
Cylinders and more...

LFY-I-5F-11

5 LPM Oxygen Concentrator

- Appealing, mobile and compact design
- Ultra-low noise



5 LPM Oxygen Concentrator







Automation is playing a pivotal role in our lives today. Likewise, automated oxygen concentrators are too out in the market to help patients better control their health. 3A Sanai Health has introduced automated 5LPM oxygen concentrators to help them get enough oxygen for the recommended time span. It, of course, varies with the respiratory condition one is suffering from.



CE 0197



Product Features



Specifications

	Nebulizer function
Name	5L Medical Oxygen Concentrator
Model	LFY-I-5F-11
Machine type	Class II, Type B
Flow rate	5L/min
Oxygen concentration	≥93% ±3%
Noise	≤50dB (A)
Voltage / Frequency	220V/50Hz, 110V/50Hz, 230V/50Hz
Power rating	400W
Output Pressure	0.03 - 0.07 Mpa
Operating Temperature	5°C - 40°C
Machine Size	520x550x200 mm
Package Size	585x635x275 mm
Net weight	21 Kg
Gross weight	22.5 Kg
Working time	24 hour uninterrupted operation
Optional configuration	Nebulizer function, remote control




SANAI HEALTH


3A Sanai Health Technology Co., Ltd.

 400 187 6068

 +86 755 82727906

 www.sanaihealth.com

 3502H, Block A, Jingji Binhe Times Square, 9289
Binhe Road, Futian District, Shenzhen, China PRC.



LOVE NATURE
LOVE SOCIETY
LOVE SOCIETY