

# PERSONAL PROTECTIVE EQUIPMENT



Face Masks, Gloves and  
Protective Clothing

# Disposable Medical Protective Clothing

Offers Complete Protection

3A Sanai has introduced a non-woven PPE suit with taping. 360-degree protection is now achievable with elastic wrists, elastic ankles, elastic hood, and the most important, the snug-fit.

It enhances breathability and offers long-lasting comfort for the medics and the paramedic's staff. Also, it allows full-coverage protection in the hazardous environment.

At the end of the day or whenever one deems fit to dispose of the PPE follow a protocol so that they don't end up touching the front surfaces of PPE, which may have been exposed to aerosols. Pulling through the sleeves to catch the inner surfaces and then rolling inside out with a pair of gloves is an ideal way to dispose of.



## Features

1. Elastic closure
2. Composed of hoods, blouses, and trousers.
3. Hydrostatic Resistant
4. Disposable
5. Does not allow self-adhesion, cracks, holes, and other defects.
6. Stitched with up to 7 to 14 stitches every 3 cm
7. Self-locking zipper
8. Easy to wear and take-off
9. Water-resistant
10. Microbiological performance

## Specifications

Available in: 160 (S), 165 (M), 170 (L), 175 (XL), 180 (XXL)

1 piece/bag, 50pcs/carton

Net Weight: 14kgs

Gross Weight: 15kgs



**SANAI HEALTH**


## **3A Sanai Health Technology Co., Ltd.**

 400 187 6068

 +86 755 82727906

 [www.sanaihealth.com](http://www.sanaihealth.com)

 3502H, Block A, Jingji Binhe Times Square, 9289  
Binhe Road, Futian District, Shenzhen, China PRC.



**LOVE NATURE**  
**LOVE SOCIETY**  
LOVE SOCIETY