

# OXYGEN SOLUTIONS



---

Oxygen Concentrators,  
Cylinders and more...

# OXYGEN CYLINDER



## Specifications

Model	Oxygen Cylinder-15	Oxygen Cylinder-20	Oxygen Cylinder-40
Volume	15 L	20 L	40 L
Outer diameter	φ159 mm	φ219 mm	φ219 mm
Height	920 mm	690 mm	1280 mm
Weight	15.6 Kg	23.9 Kg	41.0 Kg
Working pressure	15 Mpa	15 Mpa	15 Mpa
Test pressure	22.5 Mpa	22.5 Mpa	22.5 Mpa
Airtight test pressure	15 Mpa	15 Mpa	15 Mpa
Main body material	37Mn	37Mn	37Mn
Service life	20 year	20 year	20 year
Wall thickness	≥3.6 mm	≥5.0 mm	≥5.0 mm
Operating environment	-20-60 °C	-20-60 °C	-20-60 °C
Implementation standards	GB5099.3 ISO 9809-3	GB5099.3 ISO 9809-3	GB5099.3 ISO 9809-3

SP-5HB

# Oxygen Inhaler Pressure Regulator

- Comes with humidifier bottle
- Size: 765x230x63mm



Can be connected to an oxygen cylinders or an oxygen concentrator

Easy installation for up to 5 patients at a time

Individual outlet flow can be controlled individually

Avoiding cross-infection

## SPARE PARTS

- Multiple models for use with




Nasal Cannula



Humidifier bottle



HEPA filter



Non-Contact Inhaler



Oxygen mask



Silicone Nasal Cannula



Connecting Pipe



Inner Filter



Nebulizer kit



Outer filter



Oxygen Bag



Oxygen inhaler




**SANAI HEALTH**


## **3A Sanai Health Technology Co., Ltd.**

 400 187 6068

 +86 755 82727906

 [www.sanaihealth.com](http://www.sanaihealth.com)

 3502H, Block A, Jingji Binhe Times Square, 9289  
Binhe Road, Futian District, Shenzhen, China PRC.



**LOVE NATURE**  
**LOVE SOCIETY**  
LOVE SOCIETY