

# OXYGEN SOLUTIONS



---

Oxygen Concentrators,  
Cylinders and more...

OXYMATE - 7BW

# Smart Household Oxygen Concentrator

- 1-7 LPM oxygen flow rate
- Remote control operation



# OXYMATE 7BW - Smart Household Oxygen Concentrator

Oxymate 7BW offers integrated Oxygen concentration and atomization in one machine, and can be easily switched at will. Its nebulizer function further relieves respiratory symptoms hence improving circulation, decreasing anxiety levels, and lowering blood pressure. Oxymate 7BW is compact, relatively lightweight and a portable oxygen machine. It allows real-time functional monitoring just by a few clicks and offers intelligent voice output.




## Product Features



Low noise level of 40 dB



93% oxygen concentration



Oxygen is moist and not dry



1-7 LPM oxygen flow rate



Intelligent voice output



Bright and big screen display

## Specifications

<b>Name</b>	Smart Household Oxygen Concentrator
<b>Model</b>	OXYMATE 7BW
<b>Machine Type</b>	Class II
<b>Flow Rate</b>	1L/min
<b>Oxygen concentration</b>	≥93% ±3%, 7 adjustable gears
<b>Noise</b>	≤40dB (A)
<b>Voltage / Frequency</b>	220V/50Hz, 110V/50Hz, 230V/50Hz
<b>Rated power</b>	120W
<b>Output pressure</b>	0.03 - 0.07 Mpa
<b>Operating Temperature</b>	5°C - 40°C
<b>Machine Size</b>	295x252x385 mm
<b>Package Size</b>	457x320x280 mm
<b>Net Weight</b>	5.8 Kg
<b>Gross weight</b>	7.5 Kg
<b>Working time</b>	24 hour uninterrupted operation



**SANAI HEALTH**


## **3A Sanai Health Technology Co., Ltd.**

 400 187 6068

 +86 755 82727906

 [www.sanaihealth.com](http://www.sanaihealth.com)

 3502H, Block A, Jingji Binhe Times Square, 9289  
Binhe Road, Futian District, Shenzhen, China PRC.



**LOVE NATURE**  
**LOVE SOCIETY**  
LOVE SOCIETY