

# OXYGEN SOLUTIONS



---

Oxygen Concentrators,  
Cylinders and more...

OXYMOV 3

# Portable Oxygen Concentrator

Small and compact

Super long battery life




## Oxymov 3 - Portable Oxygen Concentrator

Our powerful and reliable portable oxygen concentrator is the next step in innovation for all user that are seeking to stay active all through the day and run daily errands easily with convenience and freedom. With battery backup of more than 3 hours . It uses zero discontinuous diffuse-type disper sion oxygen production to maintain a high concentration of oxygen with every breath.




### Product Features




Portable and Lightweight




High oxygen concentration



4 adjustable gears



3 hours continuous  
oxygen production



Alloy material



Large capacity battery

### Specifications

Name	Portable Oxygen Concentrator
Model	Oxymov 3
Machine type	Class I
Flow rate	1L/min
Oxygen concentration	≥93±3 %, 4 adjustable gears
Noise	≤57dB (A)
Voltage / Frequency	AC100、50/60Hz
Power rating	54W
Power Adaptation	DC 12~18V
Battery capacity	9000 mAh
Working time	Approx. 3 hours
Machine Size	178x85x227 mm
Package size	220x95x250 mm
Net weight	3 Kg
Gross weight	5 Kg



**SANAI HEALTH**


## **Shenzhen Sanai Health Technology Co., Ltd.**

 400 187 6068

 +86 755 82727906

 [www.sanaihealth.com](http://www.sanaihealth.com)

 3502H, Block A, Jingji Binhe Times Square, 9289  
Binhe Road, Futian District, Shenzhen, China PRC.



**LOVE NATURE**  
**LOVE SOCIETY**  
LOVE SOCIETY